




## Product Specification

Rack, PC, format: 12 x 4, suitable for tubes,  
S-Monovette® 15 mm and 13 mm Ø

Page 1



### Product description

Order number	93.1431
Product description	Rack, material: PC, transparent, format: 12 x 4, (LxWxH): 257 x 90 x 40 mm, for 48 collection tubes, suitable for tubes, S-Monovette® 15 mm and 13 mm Ø, 1 piece(s)/case

### Product characteristics

Suitable for	tubes, S-Monovette® 15 mm and 13 mm Ø
Storage capacity	48
Format	12 x 4
Autoclavable (max. °C)	121 °C

### Size

Width of product	90 mm
Height of product	40 mm
Length of product	257 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rack, PC, format: 12 x 4, suitable for tubes,**  
**S-Monovette® 15 mm and 13 mm Ø**

**Page 2**

**Material & colours**

Product material	Polycarbonate (PC)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	40
Piece(s) / pallet	1200
Depth of box	92 mm
Width of box	42 mm
Height of box	259 mm
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.1968 kg
Weight of case	9.61 kg
EAN of inner box	4038917086615
EAN case	4038917005159

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**