




**Product Specification**  
**Rack, PP, format: 6 x 3, suitable for 50 ml**  
**centrifuge tubes**

**Page 1**



**Product description**

Order number	93.1707
Product description	Rack, material: PP, blue, format: 6 x 3, (LxWxH): 295 x 114 x 65 mm, for 18 collection tubes, suitable for 50 ml centrifuge tubes, 1 piece(s)/bag

**Product characteristics**

Suitable for	50 ml centrifuge tubes
Storage capacity	18
Format	6 x 3
Autoclavable (max. °C)	121 °C

**Size**

Width of product	114 mm
Height of product	65 mm
Length of product	295 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rack, PP, format: 6 x 3, suitable for 50 ml**  
**centrifuge tubes**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	blue

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	540
Depth of case	125 mm
Width of case	75 mm
Height of case	300 mm
Case volume	0.0028 cbm
Weight of product	0.2298 kg
Weight of case	2.81 kg
EAN of inner box	4038917347501
EAN case	4038917347495

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**