




Product Specification

Rack, PC, format: 12 x 4, suitable for Reaction tubes 8 mm Ø e.g. REF 72.699, 72.735

Page 1



Product description

Order number	93.48.011
Product description	Rack, material: PC, transparent, format: 12 x 4, (LxWxH): 257 x 90 x 40 mm, for 48 collection tubes, suitable for Reaction tubes 8 mm Ø e.g. REF 72.699, 72.735, 1 piece(s)/case

Product characteristics

Suitable for	Reaction tubes 8 mm Ø e.g. REF 72.699, 72.735
Storage capacity	48
Format	12 x 4
Autoclavable (max. °C)	121 °C

Size

Width of product	90 mm
Height of product	40 mm
Length of product	257 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Rack, PC, format: 12 x 4, suitable for Reaction tubes 8 mm Ø e.g. REF 72.699, 72.735

Page 2

Material & colours

Product material Polycarbonate (PC)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 40
Piece(s) / pallet 1200
Depth of box 92 mm
Width of box 42 mm
Height of box 259 mm
Depth of case 398 mm
Width of case 398 mm
Height of case 340 mm
Case volume 0.0539 cbm
Weight of product 0.2426 kg
Weight of case 11.44 kg
EAN of inner box 4038917095297
EAN case 4038917007580

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT