



Product description

Order number	93.844.100
Product description	Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 257 x 62 x 55 mm, for 20 collection tubes, suitable for tubes, all S-Monovette® diameters, 1 piece(s)/case

Product characteristics

Suitable for	tubes, all S-Monovette® diameters
Storage capacity	20
Format	10 x 2
Autoclavable (max. °C)	121 °C

Size

Width of product	62 mm
Height of product	55 mm
Length of product	257 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polycarbonate (PC)
Colour of product	transparent

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	42
Piece(s) / pallet	1260
Depth of box	65 mm
Width of box	58 mm
Height of box	260 mm
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.1907 kg
Weight of case	10.14 kg
EAN of inner box	4038917121712
EAN case	4038917017961

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

