




Product Specification

Rack, PC, format: 10 x 2, suitable for tubes,
S-Monovette® 13 mm and 11 mm Ø

Page 1



Product description

Order number	93.846.100
Product description	Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 257 x 62 x 40 mm, for 20 collection tubes, suitable for tubes, S-Monovette® 13 mm and 11 mm Ø, 1 piece(s)/case

Product characteristics

Suitable for	tubes, S-Monovette® 13 mm and 11 mm Ø
Storage capacity	20
Format	10 x 2
Autoclavable (max. °C)	121 °C

Size

Width of product	62 mm
Height of product	40 mm
Length of product	257 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, PC, format: 10 x 2, suitable for tubes,
S-Monovette® 13 mm and 11 mm Ø

Page 2

Material & colours

Product material	Polycarbonate (PC)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	60
Piece(s) / pallet	1800
Depth of box	65 mm
Width of box	43 mm
Height of box	260 mm
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.17 kg
Weight of case	11.99 kg
EAN of inner box	4038917121729
EAN case	4038917017978

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT