



Product Specification

Rack, PC, format: 6 x 3, suitable for tubes,
S-Monovette® 15 mm and 13 mm Ø

Page 1



Product description

Order number	93.850
Product description	Rack, material: PC, transparent, format: 6 x 3, (LxWxH): 137 x 70 x 40 mm, for 18 collection tubes, suitable for tubes, S-Monovette® 15 mm and 13 mm Ø, 1 piece(s)/case

Product characteristics

Suitable for	tubes, S-Monovette® 15 mm and 13 mm Ø
Storage capacity	18
Format	6 x 3
Autoclavable (max. °C)	121 °C

Size

Width of product	70 mm
Height of product	40 mm
Length of product	137 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, PC, format: 6 x 3, suitable for tubes,
S-Monovette® 15 mm and 13 mm Ø

Page 2

Material & colours

Product material	Polycarbonate (PC)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	36
Piece(s) / pallet	2376
Depth of box	72 mm
Width of box	42 mm
Height of box	138 mm
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0969 kg
Weight of case	4.55 kg
EAN of inner box	4038917084369
EAN case	4038917003223

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT