



Product description

| | |
|---------------------|---|
| Order number | 93.852.300 |
| Product description | Double block rack D17, for 100 collection tubes, Ø opening: 17 mm, format: 5 x 20, grey, material: PP |

Product characteristics

| | |
|------------------------|--------|
| Number of wells | 100 |
| Storage capacity | 100 |
| Format | 5 x 20 |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|---------------------|--------|
| Diameter of opening | 17 mm |
| Width of product | 218 mm |
| Height of product | 45 mm |
| Length of product | 209 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Double block rack D17, Ø opening: 17 mm, 5 x 20,
grey**

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product grey

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of inner packaging -
Piece(s) / inner box 1
Piece(s) / case 12
Piece(s) / pallet 528
Depth of case 453 mm
Width of case 433 mm
Height of case 165 mm
Case volume 0.0324 cbm
Weight of product 0.3917 kg
Weight of case 5.41 kg
EAN of inner box 4038917122221
EAN case 4038917018296

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-05-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT