



Product description

Order number	93.854.121
Product description	Block Rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, red, material: PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	12 mm
Width of product	80 mm
Height of product	45 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	red

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	1560
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0.018 cbm
Weight of product	0.1214 kg
Weight of case	2.8 kg
EAN of inner box	4038917122641
EAN case	4038917018487

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

