



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 93.854  |
| Product description | Block Rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, grey, material: PP |

**Product characteristics**

|                        |            |
|------------------------|------------|
| Label/ Print           | with print |
| Number of wells        | 50         |
| Storage capacity       | 50         |
| Format                 | 5 x 10     |
| Autoclavable (max. °C) | 121 °C     |

**Size**

|                     |        |
|---------------------|--------|
| Diameter of opening | 12 mm  |
| Width of product    | 80 mm  |
| Height of product   | 45 mm  |
| Length of product   | 152 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | grey               |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | no                      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1             |
| Type of smallest subpackaging | -             |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 20            |
| Piece(s) / pallet             | 1560          |
| Depth of case                 | 428 mm        |
| Width of case                 | 323 mm        |
| Height of case                | 130 mm        |
| Case volume                   | 0.018 cbm     |
| Weight of product             | 0.1214 kg     |
| Weight of case                | 2.8 kg        |
| EAN of inner box              | 4038917084406 |
| EAN case                      | 4038917003247 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

