




**Product Specification**  
**Rack, PC, format: 12 x 4, suitable for reaction**  
**tubes 2 ml, Microvette®**

**Page 1**



**Product description**

Order number	93.868
Product description	Rack, material: PC, transparent, format: 12 x 4, (LxWxH): 257 x 90 x 40 mm, for 48 collection tubes, suitable for reaction tubes 2 ml, Microvette® , 1 piece(s)/case

**Product characteristics**

Suitable for	reaction tubes 2 ml, Microvette®
Storage capacity	48
Format	12 x 4
Autoclavable (max. °C)	121 °C

**Size**

Width of product	90 mm
Height of product	40 mm
Length of product	257 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rack, PC, format: 12 x 4, suitable for reaction**  
**tubes 2 ml, Microvette®**

**Page 2**

**Material & colours**

Product material	Polycarbonate (PC)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	40
Piece(s) / pallet	1200
Depth of box	92 mm
Width of box	42 mm
Height of box	259 mm
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.2411 kg
Weight of case	11.38 kg
EAN of inner box	4038917084451
EAN case	4038917003261

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-02

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**