



Product description

Order number	95.64.250
Product description	Rack, material: styrofoam, white, format: 20 x 5, (LxWxH): 416 x 107 x 50 mm, for 100 collection tubes, suitable for tubes Ø 15 and 16 mm, 1 piece(s)/blister

Product characteristics

Suitable for	tubes Ø 15 and 16 mm
Storage capacity	100
Format	20 x 5

Size

Width of product	107 mm
Height of product	50 mm
Length of product	416 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	styrofoam
Colour of product	white
Closure material	styrofoam
Colour of cap	white

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	28
Piece(s) / pallet	448
Depth of case	498 mm
Width of case	428 mm
Height of case	380 mm
Case volume	0.081 cbm
Weight of product	0.0309 kg
Weight of case	1.65 kg
EAN of inner box	4038917095587
EAN case	4038917007764

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

