



### **Product description**

Order number	95.64.251
Product description	Rack, material: styrofoam, white, format: 10 x 5, (LxWxH): 205 x 127 x 46 mm, for 50 collection tubes, suitable for 15 ml tubes, 1 piece(s)/blister

### **Product characteristics**

Suitable for	15 ml tubes
Storage capacity	50
Format	10 x 5

### **Size**

Width of product	127 mm
Height of product	46 mm
Length of product	205 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	styrofoam
Colour of product	white
Closure material	styrofoam
Colour of cap	white

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	21
Type of inner packaging	blister
Piece(s) / inner box	1
Piece(s) / case	21
Piece(s) / pallet	882
Depth of case	388 mm
Width of case	378 mm
Height of case	235 mm
Case volume	0.0345 cbm
Weight of product	0.0271 kg
Weight of case	0.58 kg
EAN of inner box	4038917095594
EAN case	4038917007771

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

