



**Product Specification**  
**CryoRack 40, PP, format: 10 x 4, suitable for**  
**CryoPure tubes**

**Page 1**



**Product description**

|                     |                                                                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 93.856.040                                                                                                                                  |
| Product description | CryoRack 40, material: PP, format: 10 x 4, (LxWxH): 100 x 210 x 22 mm, for 40 collection tubes, suitable for CryoPure tubes, 1 piece(s)/bag |

**Product characteristics**

|                  |                |
|------------------|----------------|
| Graduation       | no             |
| Suitable for     | CryoPure tubes |
| Storage capacity | 40             |
| Format           | 10 x 4         |

**Size**

|                   |        |
|-------------------|--------|
| Width of product  | 210 mm |
| Height of product | 22 mm  |
| Length of product | 100 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material                      Polypropylene (PP)

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                  no

**Packaging**

Minimum order qty.                      1  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      1  
Piece (s) / outer case                      10  
Piece(s) / pallet                            1650  
Depth of case                              228 mm  
Width of case                                218 mm  
Height of case                               155 mm  
Case volume                                0.0077 cbm  
Weight of product                        0.0679 kg  
Weight of case                              0.9 kg  
EAN of inner box                          4038917122788  
EAN case                                    4038917077293

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

