



Product description

| | |
|---------------------|--|
| Order number | 93.872.225 |
| Product description | Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, red, slip-on lid with ventilation function, cap: transparent, (LxWxH): 75 x 75 x 52 mm, format: 5 x 5, for 25 collection tubes, for CryoPure tubes 1.2?2.0 ml internal and external thread, 5 piece(s)/bag |

Product characteristics

| | |
|------------------|--|
| Graduation | no |
| Cap | slip-on lid with ventilation function |
| Closure type | slip-on lid |
| Design | with numerical coding at each aperture |
| Suitable for | CryoPure tubes 1.2?2.0 ml internal and external thread |
| Storage capacity | 25 |
| Format | 5 x 5 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|-------|
| Width of product | 75 mm |
| Height of product | 52 mm |
| Length of product | 75 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polycarbonate (PC) |
| Colour of product | red |
| Colour of cap | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 5 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 5 |
| Piece (s) / outer case | 20 |
| Piece(s) / pallet | 2640 |
| Depth of case | 156 mm |
| Width of case | 143 mm |
| Height of case | 502 mm |
| Case volume | 0.0112 cbm |
| Weight of product | 0.058 kg |
| Weight of case | 1.41 kg |
| EAN of inner box | 4038917122795 |
| EAN case | 4038917420167 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

