



Product description

Order number	93.872.625
Product description	Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, blue, slip-on lid with ventilation function, cap: transparent, (LxWxH): 75 x 75 x 52 mm, format: 5 x 5, for 25 collection tubes, for CryoPure tubes 1.2 - 2.0 ml internal and external thread, 5 piece(s)/bag

Product characteristics

Graduation	no
Cap	slip-on lid with ventilation function
Closure type	slip-on lid
Design	with numerical coding at each aperture
Suitable for	CryoPure tubes 1.2 - 2.0 ml internal and external thread
Storage capacity	25
Format	5 x 5

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Width of product	75 mm
Height of product	52 mm
Length of product	75 mm

Material & colours

Product material	Polycarbonate (PC)
Colour of product	blue
Colour of cap	transparent

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	20
Piece(s) / pallet	2640
Depth of box	262 mm
Width of box	75 mm
Height of box	75 mm
Depth of case	502 mm
Width of case	156 mm
Height of case	143 mm
Case volume	0.0112 cbm
Weight of product	0.058 kg
Weight of case	1.41 kg
EAN of inner box	4038917122832
EAN case	4038917420143

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

