



Product description

| | |
|---------------------|--|
| Order number | 93.875.281 |
| Product description | Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, red, slip-on lid with ventilation function, cap: transparent, (LxWxH): 132 x 132 x 95 mm, format: 9 x 9, for 81 collection tubes, for CryoPure tubes 3.5?5.0 ml internal and external thread, 5 piece(s)/bag |

Product characteristics

| | |
|------------------|--|
| Graduation | no |
| Cap | slip-on lid with ventilation function |
| Closure type | slip-on lid |
| Design | with numerical coding at each aperture |
| Suitable for | CryoPure tubes 3.5?5.0 ml internal and external thread |
| Storage capacity | 81 |
| Format | 9 x 9 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|--------|
| Width of product | 132 mm |
| Height of product | 95 mm |
| Length of product | 132 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polycarbonate (PC) |
| Colour of product | red |
| Colour of cap | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 5 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 5 |
| Piece (s) / outer case | 20 |
| Piece(s) / pallet | 400 |
| Depth of case | 305 mm |
| Width of case | 305 mm |
| Height of case | 508 mm |
| Case volume | 0.0472 cbm |
| Weight of product | 0.2437 kg |
| Weight of case | 5 kg |
| EAN of inner box | 4038917122948 |
| EAN case | 4038917076647 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

