



Product description

Order number	93.879
Product description	Rack, material: PP, multi-coloured, format: 12 x 8, (LxWxH): 108 x 108 x 257 mm, suitable for micro tubes and tubes with Ø 4-30 mm, 1 piece(s)/bag

Product characteristics

Suitable for	micro tubes and tubes with Ø 4-30 mm
Format	12 x 8
Autoclavable (max. °C)	121 °C

Size

Width of product	108 mm
Height of product	257 mm
Length of product	108 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, PP, format: 12 x 8, suitable for micro tubes
and tubes with Ø 4-30 mm

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	multi-coloured

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	4
Piece(s) / pallet	384
Depth of case	216 mm
Width of case	254 mm
Height of case	241 mm
Case volume	0.0132 cbm
Weight of product	0.8348 kg
Weight of case	3.6 kg
EAN of inner box	4038917355162
EAN case	4038917355155

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT