



### **Product description**

Order number	93.880
Product description	Storage box, slip-on lid, material: PP, green, format: 6 x 6, for 36 collection tubes, suitable for 15 ml tubes and tubes up to Ø 17 mm, 2 piece(s)/bag

### **Product characteristics**

Closure type	slip-on lid
Suitable for	15 ml tubes and tubes up to Ø 17 mm
Storage capacity	36
Format	6 x 6
Autoclavable (max. °C)	121 °C

### **Size**

Width of product	144 mm
Height of product	127 mm
Length of product	144 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	green
Closure material	Polypropylene (PP)
Colour of cap	transparent

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	no

### **Packaging**

Minimum order qty.	2
Type of smallest subpackaging	bag
Piece(s) / inner box	2
Piece (s) / outer case	12
Piece(s) / pallet	420
Depth of box	145 mm
Width of box	145 mm
Height of box	255 mm
Depth of case	452 mm
Width of case	251 mm
Height of case	300 mm
Case volume	0.034 cbm
Weight of product	0.2228 kg
Weight of case	3.2 kg
EAN of inner box	4038917355193
EAN case	4038917355186

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

