



Product description

Order number	95.64.922
Product description	Storage box, slip-on lid, material: cardboard, white, format: 7 x 7, for 49 collection tubes, suitable for tubes up to max. 22 mm Ø, tube height between 88.5-129 mm, 1 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Suitable for	tubes up to max. 22 mm Ø, tube height between 88.5-129 mm
Storage capacity	49
Format	7 x 7

Size

Width of product	154 mm
Height of product	130 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



Material & colours

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	27
Piece(s) / pallet	324
Depth of box	154 mm
Width of box	154 mm
Height of box	129 mm
Depth of case	550 mm
Width of case	490 mm
Height of case	415 mm
Case volume	0.1118 cbm
Weight of product	0.22 kg
Weight of case	7.3 kg
EAN of inner box	4038917095709
EAN case	4038917067447

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

