




Product Specification
**Storage box, slip-on lid, cardboard, format: 7 x 7,
for 49 collection tubes**

Page 1



Product description

Order number	95.64.994
Product description	Storage box, slip-on lid, material: cardboard, white, format: 7 x 7, for 49 collection tubes, suitable for tubes Ø 16-17 mm

Product characteristics

Closure type	slip-on lid
Suitable for	tubes Ø 16-17 mm
Storage capacity	49
Format	7 x 7

Size

Width of product	135 mm
Height of product	70 mm
Length of product	135 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	4
Piece(s) / inner box	4
Piece (s) / outer case	76
Piece(s) / pallet	960
Case volume	0.1212 cbm
Weight of product	0.1083 kg
Weight of case	8.27 kg
EAN of inner box	4038917095884
EAN case	4038917077484

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

