



Product description

Order number 65.108
Product description Push cap, white, suitable for 2 ml sample tube 73.663,
1,000 piece(s)/bag

Product characteristics

Cap flat cap
Closure type push cap
Suitable for 2 ml sample tube 73.663

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product white

Purity & certification

Product category no medical device | IVD
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 10000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 10000 |
| Piece(s) / pallet | 640000 |
| Depth of case | 398 mm |
| Width of case | 298 mm |
| Height of case | 215 mm |
| Case volume | 0.0255 cbm |
| Weight of product | 0.0004 kg |
| Weight of case | 4.1 kg |
| EAN of inner box | 4038917080866 |
| EAN case | 4038917000048 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

