



### **Product description**

Order number	65.108
Product description	Push cap, white, suitable for 2 ml sample tube 73.663, 1,000 piece(s)/bag

### **Product characteristics**

Cap	flat cap
Closure type	push cap
Suitable for	2 ml sample tube 73.663

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	white

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Push cap, white, suitable for 2 ml sample tube**  
**73.663**

**Page 2**

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	640000
Depth of box	330 mm
Width of box	295 mm
Height of box	20 mm
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0004 kg
Weight of case	4.1 kg
EAN of inner box	4038917080866
EAN case	4038917000048

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-03-10

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**