



### Product description

Order number	65.714
Product description	Screw cap, natural, suitable for tubes Ø 16.5 mm, 1,000 piece(s)/bag

### Product characteristics

Cap	tall cap
Closure type	screw cap
Suitable for	tubes Ø 16.5 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

### Purity & certification

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	265 mm
Case volume	0.042 cbm
Weight of product	0.0013 kg
Weight of case	7.14 kg
EAN of inner box	4038917083638
EAN case	4038917002530

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

