



**Product Specification**  
**Screw cap, natural, suitable for tubes and**  
**S-Monovette® Ø 13 mm**

**Page 1**



**Product description**

Order number	65.728
Product description	Screw cap, natural, suitable for tubes and S-Monovette® Ø 13 mm, 1,000 piece(s)/bag

**Product characteristics**

Cap	standard screw cap
Closure type	screw cap
Suitable for	tubes and S-Monovette® Ø 13 mm

**Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	natural

**Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap, natural, suitable for tubes and**  
**S-Monovette® Ø 13 mm**

**Page 2**

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	330000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0007 kg
Weight of case	3.76 kg
EAN of inner box	4038917083768
EAN case	4038917028288

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**