



Product Specification
Screw cap, natural, suitable for tubes and
S-Monovette® Ø 13 mm

Page 1



Product description

Order number 65.728
Product description Screw cap, natural, suitable for tubes and S-Monovette® Ø 13 mm, 1,000 piece(s)/bag

Product characteristics

Closure type screw cap
Suitable for tubes and S-Monovette® Ø 13 mm

Material & colours

Product material High Density Polyethylene (HD-PE)
Colour of product natural

Purity & certification

Product category no medical device | IVD
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Screw cap, natural, suitable for tubes and
S-Monovette® Ø 13 mm

Page 2

Packaging

Minimum order qty.	5000
Type of inner packaging	bag
Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	330000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0007 kg
Weight of case	3.76 kg
EAN of inner box	4038917083768
EAN case	4038917028288

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT