



Product description

Order number	65.729.100
Product description	Screw cap, natural, suitable for tubes Ø 15.3 mm, 100 piece(s)/bag

Product characteristics

Type of cap	screw cap
Suitable for	tubes Ø 15.3 mm

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	natural

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	5000
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	5000
Piece(s) / pallet	210000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0,0367 cbm
Weight of product	0,001 kg
Weight of case	5,62 kg
EAN of inner box	4038917119221
EAN case	4038917080101

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

