



Product Specification
**Push cap, violet, suitable for tubes Ø 15.5, 16,
16.5, 16.8 and 17 mm**

Page 1



Product description

Order number	65.793.520
Product description	Push cap, violet, suitable for tubes Ø 15.5, 16, 16.5, 16.8 and 17 mm, 1,000 piece(s)/bag

Product characteristics

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 15.5, 16, 16.5, 16.8 and 17 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	violet

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Push cap, violet, suitable for tubes Ø 15.5, 16,
16.5, 16.8 and 17 mm**

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	265 mm
Case volume	0.042 cbm
Weight of product	0.0014 kg
Weight of case	7.56 kg
EAN of inner box	4038917120906
EAN case	4038917017527

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

