



### **Product description**

Order number	65.802
Product description	Push cap, natural, suitable for tubes Ø 15.7 mm, flat, 500 piece(s)/bag

### **Product characteristics**

Cap	flat cap
Closure type	push cap
Suitable for	tubes Ø 15.7 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	10000
Piece(s) / pallet	360000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0006 kg
Weight of case	6.51 kg
EAN of inner box	4038917084185
EAN case	4038917003049

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

