



### **Product description**

Order number	65.803.005
Product description	Push cap, violet, suitable for tubes Ø 15.7 mm, 1,000 piece(s)/bag

### **Product characteristics**

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 15.7 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	violet

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0.0396 cbm
Weight of product	0.0013 kg
Weight of case	7.12 kg
EAN of inner box	4038917121088
EAN case	4038917017657

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

