



Product description

Order number	65.803.305
Product description	Push cap, violet, suitable for tubes Ø 15.7 mm, 1,000 piece(s)/bag

Product characteristics

Closure type	push cap
Suitable for	tubes Ø 15.7 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	violet

Purity & certification

Product category	no medical device IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	200000
Case volume	0.0382 cbm
Weight of product	0.0013 kg
Weight of case	7.04 kg
EAN of inner box	4038917368698
EAN case	4038917368681

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

