



### Product description

Order number	65.803
Product description	Push cap, natural, suitable for tubes Ø 15.7 mm & cuvette 67.743, 1,000 piece(s)/bag

### Product characteristics

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 15.7 mm & cuvette 67.743

### Size

Diameter	16.35 mm
Height of product	20.45 mm

### Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      5000  
Type of inner packaging                  bag  
Piece(s) / inner box                      1000  
Piece(s) / case                              5000  
Piece(s) / pallet                            210000  
Depth of case                               398 mm  
Width of case                                398 mm  
Height of case                               250 mm  
Case volume                                 0.0396 cbm  
Weight of product                         0.0013 kg  
Weight of case                              7.1 kg  
EAN of inner box                          4038917084192  
EAN case                                      4038917003056

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

