



### **Product description**

|                     |                                                                |
|---------------------|----------------------------------------------------------------|
| Order number        | 65.806.001                                                     |
| Product description | Push cap, pink, suitable for tubes Ø 13 mm, 1,000 piece(s)/bag |

### **Product characteristics**

|              |               |
|--------------|---------------|
| Closure type | push cap      |
| Suitable for | tubes Ø 13 mm |

### **Material & colours**

|                   |                                  |
|-------------------|----------------------------------|
| Product material  | Low Density Polyethylene (LD-PE) |
| Colour of product | pink                             |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 1000          |
| Piece (s) / outer case        | 10000         |
| Piece(s) / pallet             | 400000        |
| Case volume                   | 0.0382 cbm    |
| Weight of product             | 0.0002 kg     |
| Weight of case                | 7.5 kg        |
| EAN of inner box              | 4038917121170 |
| EAN case                      | 4038917033824 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

