



Product description

Order number	65.806.511
Product description	Push cap, violet, suitable for tubes Ø 13 mm, 1,000 piece(s)/bag

Product characteristics

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 13 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	violet

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	420000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0.0396 cbm
Weight of product	0.0007 kg
Weight of case	7.36 kg
EAN of inner box	4038917121255
EAN case	4038917017749

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

