



Product Specification
Push cap, natural, suitable for tubes Ø 11.5 and 12 mm

Page 1



Product description

| | |
|---------------------|--|
| Order number | 65.809 |
| Product description | Push cap, natural, suitable for tubes Ø 11.5 and 12 mm, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------|------------------------|
| Cap | standard cap |
| Closure type | push cap |
| Suitable for | tubes Ø 11.5 and 12 mm |

Size

| | |
|-------------------|---------|
| Diameter | 13.2 mm |
| Height of product | 15.1 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Low Density Polyethylene (LD-PE) |
| Colour of product | natural |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Push cap, natural, suitable for tubes Ø 11.5 and 12 mm

Page 2

Purity & certification

| | |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 10000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 10000 |
| Piece(s) / pallet | 420000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 250 mm |
| Case volume | 0.0396 cbm |
| Weight of product | 0.0006 kg |
| Weight of case | 6.66 kg |
| EAN of inner box | 4038917084239 |
| EAN case | 4038917003094 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT