



Product Specification
**Sample tube, NaCl, (LxØ): 101 x 16.5 mm, with
paper label**

Page 1



Product description

Order number	48.391
Product description	Sample tube, preparation: NaCl, screw cap, (LxØ): 101 x 16.5 mm, round base with skirted base, with paper label, 100 piece(s)/bag

Product characteristics

Type of preparation	NaCl
Preparation concentration	0.9% NaCl
Label/ Print	with paper label
Colour of print/label	black
Type of cap	screw cap
Base shape	round base with skirted base

Size

Diameter	16,5 mm
Length of product	101 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-11-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Material & colours

Product material	Polypropylene (PP)
Cap material	High Density Polyethylene (HD-PE)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	22000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0,0306 cbm
Weight of product	0,0077 kg
Weight of case	4,4 kg
EAN of inner box	4038917556613
EAN case	4038917556606

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

