



Product description

Order number	48.392.003
Product description	PP, preparation: NaCl, screw cap, (LxØ): 90 x 13 mm, round base, with paper label, 100 piece(s)/bag

Product characteristics

Type of preparation	NaCl
Preparation concentration	0.9% NaCl
Label/ Print	with paper label
Closure type	screw cap
Base shape	round base

Size

Diameter	13 mm
Length of product	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	36000
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	0.0065 kg
Weight of case	3.65 kg
EAN of inner box	4038917573870
EAN case	4038917573863

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

