



Product description

| | |
|---------------------|---|
| Order number | 55.1578 |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, suitable for flow cytometer FACSCan, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|------------------------|
| Graduation | no |
| Closure type | push cap |
| Suitable for | flow cytometer FACSCan |
| Base shape | round base |
| Volume of work | 5 ml |
| Centrifugation max (RCF) | 3000 x g |

Size

| | |
|-------------------|-------|
| Diameter | 12 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 12 mm, PP

Page 2

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 72000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 265 mm |
| Case volume | 0.042 cbm |
| Weight of product | 0.0021 kg |
| Weight of case | 4.71 kg |
| EAN of inner box | 4038917087087 |
| EAN case | 4038917005555 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT