



Product Specification
Tube, 5 ml, (LxØ): 75 x 12 mm, PS

Page 1



Product description

Order number	55.1579
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, suitable for flow cytometer like FACSCanto, FACS Aria, FACSCalibur, LRS2, FC500, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	flow cytometer like FACSCanto, FACS Aria, FACSCalibur, LRS2, FC500
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 12 mm, PS

Page 2

Size

Diameter	12 mm
Length of product	75 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	112000
Depth of box	394 mm
Width of box	77 mm
Height of box	179 mm
Depth of case	428 mm
Width of case	383 mm
Height of case	175 mm
Case volume	0.0287 cbm
Weight of product	0.0026 kg
Weight of case	5.82 kg
EAN of inner box	4038917157735
EAN case	4038917157728

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT