



### Product description

Order number	55.460
Product description	Tube, working volume: 16.5 ml, (LxØ): 130 x 17 mm, material: PS, round base, transparent, push cap, 250 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	16,5 ml
Centrifugation max (RCF)	3500 x g

### Size

Diameter	17 mm
Length of product	130 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Tube, 16.5 ml, (LxØ): 130 x 17 mm, PS**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0086 kg
Weight of case	9.35 kg
EAN of inner box	4038917081689
EAN case	4038917000796

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**