



Product description

Order number	55.461
Product description	Tube, working volume: 14 ml, (LxØ): 119 x 15.5 mm, material: PS, round base, transparent, 500 piece(s)/bag

Product characteristics

Graduation	no
Type of cap	push cap
Base shape	round base
Volume of work	14 ml
Centrifugation max (RCF)	3500 x g

Size

Diameter	15.5 mm
Length of product	119 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of inner packaging bag
Piece(s) / inner box 500
Piece(s) / case 1000
Piece(s) / pallet 30000
Case volume 0.0539 cbm
Weight of product 0.0058 kg
Weight of case 6.29 kg
EAN of inner box 4038917081696
EAN case 4038917000802

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

