



### Product description

Order number	55.463
Product description	Tube, working volume: 14 ml, (LxØ): 105 x 16.8 mm, material: PS, round base, transparent, push cap, 250 piece(s)/StackPack

### Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	14 ml
Centrifugation max (RCF)	3500 x g

### Size

Diameter	16.8 mm
Length of product	105 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of box	360 mm
Width of box	107 mm
Height of box	184 mm
Depth of case	388 mm
Width of case	378 mm
Height of case	235 mm
Case volume	0.0345 cbm
Weight of product	0.0062 kg
Weight of case	6.95 kg
EAN of inner box	4038917081719
EAN case	4038917000826

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

