



### **Product description**

Order number	55.466.005
Product description	Tube, working volume: 11.5 ml, (LxØ): 100 x 15.7 mm, material: PS, round base, transparent, push cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	11,5 ml
Centrifugation max (RCF)	3500 x g

### **Size**

Diameter	15.7 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0055 kg
Weight of case	6.19 kg
EAN of inner box	4038917106795
EAN case	4038917011563

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

