



### **Product description**

Order number	55.466
Product description	Tube, working volume: 11.5 ml, (LxØ): 100 x 15.7 mm, material: PS, round base, transparent, push cap, 250 piece(s)/StackPack

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	11,5 ml
Centrifugation max (RCF)	3500 x g

### **Size**

Diameter	15.7 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of box	336 mm
Width of box	102 mm
Height of box	174 mm
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0055 kg
Weight of case	6.1 kg
EAN of inner box	4038917081740
EAN case	4038917000857

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

