



Product Specification
Tube, 26 ml, (LxØ): 87 x 23.5 mm, PS

Page 1



Product description

Order number	55.471.004
Product description	Tube, working volume: 26 ml, (LxØ): 87 x 23.5 mm, material: PS, round base, transparent, push cap, natural, cap assembled, sterile, 4 piece(s)/bag

Product characteristics

Graduation	no
Cap	standard cap
Closure type	push cap
Base shape	round base
Volume of work	26 ml
Centrifugation max (RCF)	3500 x g

Size

Diameter	23.5 mm
Length of product	87 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	400
Type of smallest subpackaging	bag
Piece(s) / inner box	4
Piece (s) / outer case	400
Piece(s) / pallet	12800
Depth of case	648 mm
Width of case	420 mm
Height of case	145 mm
Case volume	0.0395 cbm
Weight of product	0.0091 kg
Weight of case	4.61 kg
EAN of inner box	4038917107211
EAN case	4038917011709

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

