



Product description

Order number	55.475
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PS, round base, transparent, push cap, without cap, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	112000
Depth of box	416 mm
Width of box	77 mm
Height of box	193 mm
Depth of case	428 mm
Width of case	383 mm
Height of case	175 mm
Case volume	0.0287 cbm
Weight of product	0.0031 kg
Weight of case	6.89 kg
EAN of inner box	4038917081832
EAN case	4038917000949

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

