



### **Product description**

Order number	55.484.001
Product description	Tube, working volume: 3.5 ml, (LxØ): 55 x 12 mm, material: PS, round base, transparent, push cap, natural, cap assembled, sterile, 100 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard cap
Closure type	push cap
Base shape	round base
Volume of work	3,5 ml
Centrifugation max (RCF)	5000 x g

### **Size**

Diameter	12 mm
Length of product	55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	natural
Tube material	Polystyrene (PS)

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	64000
Depth of case	648 mm
Width of case	420 mm
Height of case	145 mm
Case volume	0.0395 cbm
Weight of product	0.0026 kg
Weight of case	5.91 kg
EAN of inner box	4038917107846
EAN case	4038917012010

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

