



Product description

Order number 55.484.005
Product description Tube, working volume: 3.5 ml, (LxØ): 55 x 12 mm, material: PS, round base, transparent, 1,000 piece(s)/bag

Product characteristics

Graduation no
Closure type push cap
Base shape round base
Volume of work 3,5 ml
Centrifugation max (RCF) 5000 x g

Size

Diameter 12 mm
Length of product 55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 4000
Type of smallest
subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 4000
Piece(s) / pallet 120000
Depth of case 398 mm
Width of case 398 mm
Height of case 340 mm
Case volume 0.0539 cbm
Weight of product 0.0019 kg
Weight of case 8.48 kg
EAN of inner box 4038917107853
EAN case 4038917012027

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

