



### **Product description**

Order number 55.484.005  
Product description Tube, working volume: 3.5 ml, (LxØ): 55 x 12 mm, material: PS, round base, transparent, push cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation no  
Closure type push cap  
Base shape round base  
Volume of work 3,5 ml  
Centrifugation max (RCF) 5000 x g

### **Size**

Diameter 12 mm  
Length of product 55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	4000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	4000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0019 kg
Weight of case	8.48 kg
EAN of inner box	4038917107853
EAN case	4038917012027

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

