



Product description

Order number	55.484
Product description	Tube, working volume: 3.5 ml, (LxØ): 55 x 12 mm, material: PS, round base, transparent, push cap, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	3,5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	12 mm
Length of product	55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 4000
Type of smallest StackPack
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 4000
Piece(s) / pallet 168000
Depth of box 381 mm
Width of box 57 mm
Height of box 176 mm
Depth of case 388 mm
Width of case 378 mm
Height of case 250 mm
Case volume 0.0367 cbm
Weight of product 0.0019 kg
Weight of case 8.66 kg
EAN of inner box 4038917081917
EAN case 4038917001021

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

