



### Product description

Order number	55.511
Product description	Tube, working volume: 14 ml, (LxØ): 119 x 15.5 mm, material: PP, round base, transparent, 500 piece(s)/bag

### Product characteristics

Graduation	no
Type of cap	push cap
Base shape	round base
Volume of work	14 ml
Centrifugation max (RCF)	3500 x g

### Size

Diameter	15.5 mm
Length of product	119 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min.
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1000
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	30000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0048 kg
Weight of case	5.43 kg
EAN of inner box	4038917082068
EAN case	4038917001151

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

